


How to Use Your Ellipta

1. Remove foil package from the box.
2. Open foil package and remove the Ellipta device.



3. Hold device upright.
4. Slide the cover open until you hear a “click”.
5. **The medication is loaded and now ready to take.**



6. Tilt your head back slightly and sit up straight.
7. Blow out and away from the device.



8. Seal your lips around the mouthpiece.
9. **Be sure your fingers do not cover the vent holes on the front of the device.**
10. Breathe in slowly & deeply.
11. Hold your breath for up to 10 seconds.



12. Wipe mouthpiece with a dry cloth.
13. Slide the cover closed until you hear a “click”.
14. When you are done, rinse, gargle and spit.



15. There is a counter on the front of the device that will count down by one each time you use it. It will read “0” when the device is empty.
16. Store in a cool, dry place.